

Dotcount User Guide

Joacim Thomassen

6th September, 2012

Chapter 1

System requirements

Dotcount was built and tested on Ubuntu 10.04 LTS. A workstation with this system up and running is thus necessary to use Dotcount as it is. If you don't already have an Ubuntu workstation you're advised to seek information and guidance from Ubuntu's own website ubuntu.com.

Dotcount depends on another software component, the NetPBM library. Ubuntu will install this component for you if you search for `libnetpbm10` in the Ubuntu Software Center and select install.

With Ubuntu and the NetPBM component installed your final step is to install Dotcount itself. This is done by copying the executable file called `dotcount` to your own system. Copy this single file to your system's local program path¹. See listing 1.1 for the command line you must enter in the Terminal program. That's it. You now have a system with Dotcount installed.

Listing 1.1: Install program with copy command. Replace MYMEDIA with the name given to your media containing dotcount.

```
sudo cp /media/MYMEDIA/bin/dotcount /usr/local/bin/
```

¹local program path = `/usr/local/bin/`

Chapter 2

Overview




Figure 2.1: Program overview

Dotcount is an image analysis program to count dots and approximate each dot's diameter. Dotcount has a command line interface and expects to be started together with the original image, a prepared version of the original image, and two scale sizes. The program then locates all the dots, marks them in a copy of the original image and reports the total amount and approximated dot diameter. This report is stored in two text files and can be used in further statistical analysis of the results. For an overview see figure 2.1.

Chapter 3

Preparing program input

Dotcount has so far been used to analyze SEM images showing quantum dots in the layer structure of a new material. The laboratory image must have meta data about the scale to give meaningful approximations about dot diameters. Dotcount do need some manual help to reveal and separate out the dots from the background and image noise. This is done by preparing the lab image for Dotcount in an standard image editor like The Gimp.




Figure 3.1: Program input preparation

Thresholds and curves are used to separate tone values in the greyscale image and must result in a binary image with only white and black values. This preparation is documented in the instructional video “Input preparation”.

In figure 3.1 an overview of the preparation procedure is illustrated. The details of this image editing routine is presented in the instructional video and are outlined in listing 3.

1. Image the surface of your sample material. Save the image as a digital image like TIFF or similar.
2. Open the lab image in an image editor like The Gimp. Do the following.
 - (a) Extract scale information.

- (b) Crop image so only the surface part that you wish to analyze is visible. No borders or meta data like logo etc. can be present.
 - (c) Increase contrast if necessary.
 - (d) Convert the image to PNM by saving it in this format. This saved image is called the original image.
3. Work on the original image (cropped PNM) to be sure you have the identical image data for the prepared image we're about to make.
- (a) Adjust levels and curves to get the cleanest dot representation.
 - (b) Save the image with a different name, like sample01prepared.pnm. This saved image is called the prepared image.

Chapter 4

Program usage

Dotcount is executed from the command line and will save its output to the same directory as it is invoked from. If you have your original image and prepared image in a directory named *sample1* you could invoke dotcount from this directory and the result will be saved in this directory. The marked image, plot file and log file can then be found in your directory *sample1* together with the relevant input data. You can also have multiple input data in the same directory as the output from dotcount is named using the name of the input file. Dotcount expects the arguments shown in listing 4.1. An example run is illustrated in listing 4.2. The result from this run is the three outputted files shown in listing 4.3.

Listing 4.1: Dotcount usage

```
dotcount <prepared.pnm> <original.pnm> <nm.range> <px.range>
```

Listing 4.2: Dotcount example

```
dotcount my_prepared.pnm my_original.pnm 100 43
```

Listing 4.3: Output from listing 4.2 example

```
my_original.marked.pnm my_original.plot my_original.log
```


Chapter 5

Program result

Running Dotcount results in three output units. The copy of the original image with red marking of dot perimeters are a good feedback on how well you prepared the image data before processing. The log file contains general information as the example in listing 5.1.

Listing 5.1: Reported logfile

```
Processed by: Dotcount 2.0 (2010-3-17)
Sample: testimage_original
mm range: 100
px range: 43
Density per px2: 0.0024846 dots/px2
Density per cm2: 4.59402e+10 dots/cm2
Dots: 1697
Total area in px2: 683008 px2
Image height: 1.55116 um2
Image width: 2.3814 um2
Total area in um2: 3.69393 um2
Average circumference in px: 26 px (~ Avg.D: 8.27606 px)
Average circumference in nm: 60.4651 nm (~ Avg.D: 19.2466 nm)
```